

Analysis of water content distribution in food by IRT-7000S

<Introduction>

Even though the detector of IRT-7000S is single-element detector, it enables high speed imaging measurement with 1/10 measurement time due to high speed mapping system as compared with the conventional IR microscope.

Table 1 shows the time required for measurement of $200 \times 200 \mu\text{m}$ area by multi-channel detector (16×1 linear array detector) and single-element detector (IRT-7000), and Figure 1 shows the image of measurement protocol by each detector. Generally multi-channel detector is used for high speed imaging measurement, and its measurement area size is determined by the magnification of Cassegrain mirror. On the other hand, the measurement area size by single-element detector is determined by the aperture size.

Accordingly, single-element detector is considered to have an advantage because the aperture size is changeable, so that for large area measurement with large aperture size, the measurement time may be shorter than multi-channel detector. Therefore, for the measurement of $200 \times 200 \mu\text{m}$ area, the measurement time by single-element detector with $50 \times 50 \mu\text{m}$ aperture size can be shorter than linear array detector with $x 32$ Cassegrain, even the measurement points are reduced.

One of the significances to measure large area in high speed is the requirement for unstable sample such as food or tissue which may denature in a short time. In this experiment, we would like to show visually the difference of water content distribution in 2 types of noodles cooked by professional and amateur. By using IRT-7000S *1, high speed imaging measurement was done for the samples in a few millimeter block before drying.

Detector	Linear array (16×1)	Single-element (IRT-7000S)	Single-element (IRT-7000S)
Aperture size	$6.25 \times 6.25 \mu\text{m}$	$25 \times 25 \mu\text{m}$	$50 \times 50 \mu\text{m}$
Measurement points	32×32	8×8	4×4
Measurement time *3	Approx. 7 sec	Approx. 7 sec	Approx. 2 sec

Table 1 Comparison of detectors and aperture size in $200 \times 200 \mu\text{m}$ area measurement *2

*1: IRT -7000S can be upgraded to IRT-7000

*2 Measurement condition: $x 32$ Cassegrain, resolution: 16 cm^{-1} , accumulation: 1

*3 Measurement time depends on measurement points




Figure 1 Comparison of measurement area in $200 \times 200 \mu\text{m}$ area

In order to compare water content in noodle cross-section, the sample was prepared quickly and measured with high speed. The noodles cooked by professional and amateur were cut in a few millimeters thickness and sandwiched by CaFs windows in order to avoid drying (Figure 2).

<Measurement condition>

Instrument: IRT-7000S
 Measurement points: 38×54 points (2052 points) (professional)
 Measurement method: transmittance
 Measurement size: 3.7 mm \times 5.3 mm (professional), 4.3 mm \times 4.6 mm (amateur)
 Aperture: 100 \times 100 μm
 Detector: MCT
 Measurement time: approx. 5 minutes (professional)
 Resolution: 16 cm^{-1}




Figure 2. Process of measurement

<Results>

The water content distribution in both noodle cooked by professional and amateur was compared by calculating the peak area ratio of starch ($3872 - 4165 \text{ cm}^{-1}$) and water ($1882 - 2321 \text{ cm}^{-1}$) (Figure 3).

Figure 3 shows that water content in the noodle cooked by amateur was rich in the center, while water content in the noodle cooked by professional was distributed in the layer between outer side and core part. In addition, as a results of sensory test, almost all testers of 20 people said that the noodle cooked by professional was more delicious. In this experiment, By high speed mapping of few millimeters sample with $100 \mu\text{m}^2$ measurement spot size, it was confirmed that the taste was related to the water content distribution,.

By using single-element detector system with changeable aperture size when combined with high speed mapping, it enables to measure the samples with area as large as few milimeters in a shot time and visualize the sample in molecular level (such as water content distribution in food and medicine).

Using such system will be very useful for the study of food/medicine, if the relationship between the obtained image and food taste or effectiveness of medicine could be determined.




Figure 3. Water content distribution of fried noodle cooked by professional and amateur