

## High-Speed Separation of ATP and its Degradation Products by Ultra High-performance Liquid Chromatography and its Application to Evaluation of Degree of Freshness of Fish Meat

### Introduction

ATP(adenosine triphosphate) in fish meat is decomposed as time elapses following the route shown below.

ATP → ADP(adenosine diphosphate) → AMP(adenosine monophosphate) → IMP(inosinic acid) →  
Ino(inosine) → Hypo(hypoxanthine)

There is a K value which indicates the degree of freshness of fish meat, and is defined as follows.

$$K \text{ value}(\%) = \frac{[\text{Ino} + \text{Hypo}]}{[\text{ATP} + \text{ADP} + \text{AMP} + \text{IMP} + \text{Ino} + \text{Hypo}]} \times 100$$

It is known that fish meat can be used for sashimi if K value is less than 20% and can be used for cooking and processing if K value is 20 - 60%.

In this report, the degree of freshness of fish meat is measured by calculating K values using Ultra High-performance Liquid Chromatography (UHPLC). In addition, the time course of K values, i.e., degradation of the freshness, was measured. K value was calculated using the free calculation function in the ChromNAV CDS.

**Keyword :** UHPLC, ATP, ADP, AMP, IMP, Ino, Hypo, 2.0 μm, C18 column, UV-Vis detector

### Experimental

#### Equipment.


Pump: X-LC 3185PU x 2  
Degasser: X-LC 3080DG  
Mixer: X-LC 3180MX  
Column oven: X-LC 3067CO  
Autosampler: X-LC 3159AS  
Detector: X-LC 3070UV

#### Conditions.

Column: X-PressPak AQ-C18-W (3.0 mmID x 50 mmL, 2.0 μm)  
Eluent A: 100 mM Phosphate buffer (pH 4.2)  
Eluent B: 100 mM Phosphate buffer (pH 4.2)/Acetonitrile (50/50)  
Gradient condition: (A/B), 0 min(100/0) → 2.2 min(100/0) → 6.0 min(50/50) → 7.0 min(50/50) → 7.05 min(100/0) 1 cycle; 10 min  
Flow rate: 0.6 mL/min  
Column temp.: 30°C  
Wavelength: 260 nm  
Injection volume: 1 μL  
Standard sample: ATP and degradation products 1.0 μg/mL each

### Results


Fig. 1 shows chromatogram of adenosine nucleotides. Analysis time was shortened approximately 1/8 times as compared with conventional HPLC without sacrificing the resolution of each component.



**Fig. 1.** Chromatogram of standard mixture of adenosine nucleotides 1: ATP, 2: ADP, 3: IMP, 4: Hypo, 5: AMP, 6: Ino

copyright©JASCO Corporation


Fig. 2 shows a chromatogram of adenosine nucleotides in sashimi grade tuna fish meat two days after purchase. Target six components are detected without any interference by contaminants. This sample's K value was 8 %.



**Fig. 2.** Chromatogram of adenosine nucleotides in sashimi grade tuna fish meat (stored in a refrigerator at 4°C for two days after purchase). 1: ATP, 2: ADP, 3: IMP, 4: Hypo, 5: AMP, 6: Ino

**Preparation.** 0.4 M perchloric acid aqueous solution(20 mL) was added to tuna fish meat (2.5 g) and homogenized. Then 2 M potassium carbonate aqueous solution(1 mL) was added to the obtained supernatant (5 mL) and applied to centrifugal separation. Obtained supernatant was filtered with 0.2 μm membrane filter.

Fig. 3 shows the time course of the degradation of tuna fish meat.



**Fig. 3.** Time course of the degradation of tuna fish meat