

## Analysis of Vitamin K using Reduction Catalyst Column

Vitamin K is essential for the synthesis of several proteins involved in blood clotting and calcium-binding. The vitamin K-dependent proteins include the bone proteins such as osteocalcin which has a function of bone metabolism, and coagulation/anticoagulation proteins which control blood coagulation. However, its physiological and physiological roles and behavior in vivo are still unclear. To clarify these, a simple and highly sensitive measurement method is needed.

As of today, the electrochemical reduction method, post-column derivatization method by sodium tetraborate solution and reduction method by platinum column are used for highly sensitive analysis of vitamin K, using reduction of the K-quinone to the fluorescent K-hydroquinone.

In this report, vitamin K, reduced by a reduction column set behind a separation column, is detected by the fluorescence detector with high selectivity and sensitivity. Fig. 1 shows the system configuration, and Fig. 2 shows the measurement principle. Chromatogram of standard samples are shown in Fig. 3.

### Conditions:

Column:	CrestPak C18S (4.6 mm I.D. x 150 mm)
Eluent:	CH <sub>3</sub> OH/C <sub>2</sub> H <sub>5</sub> OH (70/30)
Flow rate:	1.0 mL/min
Column temperature:	40 degree celsius
Reduction catalyst column:	CATALYSISPAK PT (4.6 mm I.D. x 10 mm)
Reaction temperature:	40 degree celsius
Wavelength:	Ex 320 nm, Em 430 nm, Gain x100
Sample:	STD mixture (Vitamin K <sub>1</sub> , K <sub>2</sub> , K <sub>3</sub> ) (1.0 mg/L each)
Injection volume:	10 μl

*Keywords:* 1. Vitamin K<sub>1</sub>, 2. Vitamin K<sub>2</sub>, 3. Vitamin K<sub>3</sub>,  
4. CrestPak C18S, 5. CATALYSISPAK PT




Fig. 1 Flow diagram




Fig. 2 Measurement principle




Fig. 3 Chromatogram of vitamin Ks