

## Analysis of indometacin using direct injection in of rat plasma

When analyzing for drugs in biological samples such as blood plasma, preparation by deprotenization is usually required.

In the present report, indomethacin in rat plasma was analyzed by online deprotenization using a deprotenization column with column switching.

Fig. 1 shows the chromatogram of rat plasma spiked with indomethacin standard (1.0 mg / ml).

**Keywords:** 1. Indometacin, 2. rat plasma, 3. CrestPakC18T-5, 4. UV, 5. column switching

### Conditions:

Pretreatment Column: CAPCELL PAK MF Ph-1,  
(4.6mm I.D. x 50mmL)  
Pretreatment Eluent: 100mM Na<sub>2</sub>HPO<sub>4</sub> (pH6.9) /  
CH<sub>3</sub>CN (87/13)  
Flow rate: 1.0mL/min  
Switching time: 5.7~7.3min  
Separation Column: CrestPak C18T-5  
(4.6mm I.D. x 250mmL)  
Eluent: 100mM Na<sub>2</sub>HPO<sub>4</sub> +20mM SDS  
(pH2.3)/CH<sub>3</sub>CN (50/50)  
Wavelength: 254nm  
Range: 0.001AU/10mV  
Column temperature: 40 degree celsius  
Injection volume: 20uL  
Sample: Indometacin STD + rat plasma  
(1.0ug/mL)





Fig. 1 chromatogram of indomethacin (standard added at 1.0 mg / ml) in rat plasma.