

Separation of Basic Drugs by Supercritical Fluid Chromatography

Introduction

It is well known that basically, SFC with a column of high polarity material shows the same retention action as normal phase chromatography, and therefore it is believed to be difficult to separate aqueous high polarity components. However, if a little volatile acid, base, or salt is added to modifier solvent (alcohol or etc.), the shape of polar component's peak can be improved, and components with long retention time can be eluted with appropriate retention time.

In this experiment, basic drugs were separated using 2-Ethylpyridine as column, Ammonium acetate in Methanol as modifier solvent.

Keyword: SFC, Basic Drugs, 5 μ L, 2-Ethylpyridine column, PDA Detector

Experimental

Equipment

CO₂ Delivery Pump: PU-2080-CO2
 Modifier Pump: PU-2080
 Mixer: MX-2080-32
 Autosampler: AS-2059-SF
 Column Oven: CO-2060
 PDA Detector: MD-2018
 (High pressure cell)
 Back Pressure Regulator: BP-2080

Conditions

Column: 2-Ethylpyridine 60A (4.6 mmID x 250 mmL, 5 μ m)
 (Princeton Chromatography Inc.)
 CO₂ Flow rate: 3.0 mL/min
 Modifier: 20 mM Ammonium acetate in Methanol
 Modifier gradient: 0 min (0.2 mL/min) \rightarrow 6 min (0.2 mL/min) \rightarrow
 13 min (1.0 mL/min) \rightarrow 18 min (1.0 mL/min) \rightarrow
 18.05 min (0.2 mL/min) 1 cycle: 30 min
 Column temp.: 40°C
 Back Pressure: 15 MPa
 Wavelength: 200-400 nm
 Injection volume: 5 μ L
 Standard sample: Mixture 0.1 mg/mL each

Result

Chromatogram and contour plot of standard mixture of basic drugs (wavelength: at 220 nm) are shown in Fig.1. Polar components such as Berberine and Maleic were also eluted.




Fig. 1. Chromatogram of standard mixture of basic drugs

1: Caffeine, 2: Hexobarbital, 3: Amitriptyline, 4: Chlorpheniramine, 5: Imipramine
 6: Quinine, 7: Atropine, 8: Acetaminophen, 9: Berberine, 10: Maleic acid

copyright©JASCO Corporation