

Analysis of Formaldehyde in Air using Post-column Derivatization by High Performance Liquid Chromatography

Introduction

It is becoming a big concern that aldehydes such as the formaldehyde and acetaldehyde as an environmental pollutant may contaminate environment such as the atmosphere, lakes and marshes, reservoirs, and rivers. Therefore, it is an object of various regulations like Air Pollution Control Law, Water Supply Law, and Offensive Odor Control Law, etc. in this country. As a method to measure aldehydes in atmosphere using HPLC, general method employed is collecting aldehydes in atmosphere by using the sampler enclosing silicagel to contain 2,4- DNPH and extracting aldehydes derivatized after collection by the acetonitrile. JASCO has introduced so far the analysis of formaldehyde and acetaldehyde by the post column fluorescence derivatization method using 1,3-Cyclohexanedione as a labeling reagent.

Here, it was examined to use the same method to analyze aldehydes in atmosphere. As a result, aldehydes in atmosphere extracted by water were analyzed successfully as well as the case of aldehydes in water.

Keyword : Aldehyde genera, Passive gas tube, 1,3-Cyclohexanedione, Post column derivatization method, Shodex RSpak KC-811 6E, Fluorescence detector

Experimental

Equipment


Eluent pump: PU-2080
 Reagent pump: PU-2085
 Degasser: DG-2080-53
 Autosampler: AS-2057
 Column oven: CO-2060
 Reaction oven: RO-2061
 Detector: FP-2020

Conditions

Column: Shodex RSpak KC-811 6E (6.0 mmID x 250 mmL)
 Eluent: 3 mM Perchloric acid
 Flow rate: 1.0 mL/min
 Reagent: 1,3-Cyclohexanedione in ammonium acetate buffer
 Reagent flow rate: 0.4 mL/min
 Column temp.: 60°C
 Reaction temp.: 120°C
 Wavelength: Ex. 366 nm, Em. 440 nm, Gain x10
 Injection volume: 50 µL
 Standard sample: Formaldehyde 0.1 mg/mL

Result

In Fig. 1, 1,3-Cyclohexanedione reaction formula of post column derivatization method is shown and Fig. 2 illustrates Flow system diagram.



Decahydroacridine-1,8-dione

Fig. 1. 1,3-Cyclohexanedione reaction formula of post column derivatization method




Fig. 2. Flow system diagram

- No.1 : Eluent
- No.2 : Degasser (DG-2080-53)
- No.3 : Pump for eluent (PU-2080)
- No.4 : Autosampler (AS-2057)
- No.5 : Column oven (CO-2060)
- No.6 : Column (Shodex RSpak KC-811 6E)
- No.7 : Reagent
- No.8 : Pump for reagent (PU-2085)
- No.9 : Reaction oven (RO-2061)
- No.10: Reaction coil
- No.11: Cooling coil
- No.12: Fluorescence detector (FP-2020)
- No.13: Backpressure coil

Fig. 3 explains the usage of passive gas tube for collecting the sample in atmosphere and Fig. 4 shows the extraction method from passive gas tube.




Fig. 3. Usage of passive gas tube for collecting the sample in atmosphere **Fig. 4.** Extraction method from passive gas tube

Fig. 5 shows the chromatogram of the sample in the laboratory atmosphere collected by passive gas tube. The calculated concentration of Formaldehyde obtained was 0.0054 mg/L.




Fig. 5. Chromatogram of sample in the laboratory atmosphere
1: Formaldehyde (0.0054 mg/L)