

## ULTRASONIC POWER DETERMINATION

Although other methods have been used to determine the amount of power that is delivered to a sample, such as using a hydrophone immersed below the probe, the most common and easiest method is outlined below.

- 1) Turn on the equipment
- 2) Set the amplitude as required
- 3) With the probe in air, not immersed in a sample, record the amount of watts displayed on the power monitor
- 4) Without changing the amplitude setting, immerse the probe into the sample and record the amount of watts displayed on the power monitor
- 5) The difference in power readings between step 3 and 4 above, is the amount of power being delivered to the sample in watts
- 6) To obtain the Power Density in watts/cm<sup>2</sup>, divide the number of watts obtained in step 5 by the area of the probe tip.

$$\text{Area} = (\text{diameter}/2)^2 \times \pi \text{ or } \pi r^2$$

See below for an example using a 3mm probe.

diameter of 3mm probe = 3mm/10 = .3cm = .15cm radius

.15<sup>2</sup>cm X 3.142 (π) = .0707cm<sup>2</sup> area of probe

Divide the # of watts delivered into the sample by the area of the probe

At 20% amplitude, with 1 watt delivered into the sample, the power density would be  $1/.0707 = 14 \text{ watt/cm}^2$

Amplitude	Output watts in Air	Output watts in Sample	Output watts Delivered to Sample	Watts/cm <sup>2</sup>
20	1	2	1	14
40	3	6	3	42
60	5	10	5	71
80	8	15	7	99

Intensity is expressed as power per surface area of an ultrasonic probe  
(watts/cm<sup>2</sup>)

**NOTE:**

The greater the resistance to the movement of the probe, the greater the amount of power that will be delivered to the probe.